



is valid. An example of a standard series and recovery of an endotoxin spike of 0.5 EU/mL in sample products is presented below.

REPRESENTATIVE LINEAR ANALYSIS

SPECIMEN	CSE (EU/mL)	MEAN ONSET TIME (sec)	RECOVERY (EU/mL)
STD 1	50.0	337	
STD 2	5.0	493	
STD 3	0.5	801	
STD 4	0.05	1400	
STD 5	0.005	3025	
AVG. SPL1	----	****	<0.005
AVG. PPC	0.5	822	0.706 (141%)

Fitted regression equation: Log(y) = -0.236\*Log(x) + 2.879  
Correlation Coefficient: r = -0.991

In this example, the positive product control (PPC) for each sample (SPL) yielded endotoxin recovery consistent with absence of inhibition. The negative control should be significantly lower than the lowest standard endotoxin concentration.

If the absolute value of the correlation coefficient is greater than 0.980, a polynomial regression model can be used to describe the standard curve. The polynomial regression is calculated for the user using Endoscan-V software<sup>12</sup>, or equivalent software. For additional information, please refer to the related software manuals and the section describing polynomial regression below.

POLYNOMIAL REGRESSION

A polynomial regression model may be fitted to the standard curve, providing the absolute value of the linear correlation coefficient is greater than, or equal to 0.980. Software packages such as Endoscan-V<sup>12</sup> are capable of fitting the best polynomial function to the data.

NOTE: A polynomial standard curve can not be used for initial qualification assays. The FDA accepts linear regression models for such assays, per the regulatory guidelines.

The order of the polynomial regression is determined as follows:

Endoscan-V Software n-1

Where n = the number of standards utilized.

The example assay (above), analyzed using a fourth degree polynomial (five standards), would yield the following results:

REPRESENTATIVE POLYNOMIAL ANALYSIS

SPECIMEN	CSE (EU/mL)	MEAN ONSET TIME (sec)	RECOVERY (EU/mL)
STD 1	50.0	337	
STD 2	5.0	493	
STD 3	0.5	801	
STD 4	0.05	1400	
STD 5	0.005	3025	
AVG. SPL1	----	****	<0.005
AVG. PPC	0.5	822	0.445 (89%)

Fitted regression equation:  
Log(y) = 0.0031\*Log(x)<sup>4</sup> + 0.0001\*Log(x)<sup>3</sup> + 0.011\*Log(x)<sup>2</sup> -0.216\*Log(x) + 2.837  
  
Linear Correlation Coefficient: r = -0.991

LIMITATIONS OF PROCEDURE

Samples may be tested by LAL methods provided that no inhibition or enhancement conditions are present that can not be eliminated by an acceptable dilution (refer to MVD calculation) or sample-pretreatment, such as buffering. If the LAL method cannot be validated at a concentration within the maximum valid dilution, the LAL test cannot be substituted for the USP Pyrogen Test.

PERFORMANCE CHARACTERISTICS

**Linearity:** The linearity of the standard curve within the concentration range used to determine endotoxin levels must be verified. No less than 3 endotoxin standards, spanning the desired concentration range, should be assayed at least in triplicate.<sup>10</sup> The absolute value of the coefficient of correlation, r, shall be greater than or equal to 0.980.<sup>10</sup>

QUALITY CONTROL PROCEDURES FOR THE KINETIC-COLORIMETRIC METHOD

Follow the U.S. Pharmacopeia, 25<sup>th</sup> revision for end-product testing using Kinetic Colorimetric methods, including the attainment of a Positive Product Control within 50 – 200% of theoretical concentration. Standard curves must have a correlation coefficient ≤ -0.980.

When consistency in standard curves is observed, an archived curve may be generated, using 3 assays, as described by the Interim Guidance.<sup>10</sup> If an archived curve is used, the positive water control must recover a mean endotoxin concentration within ± 25% of this point on the archived curve.<sup>10</sup>

Charles River Endosafe has developed a guide for initial qualification of kinetic incubating microplate readers.

EXPECTED VALUES

Endochrome-K™ LAL Reagent is standardized against the U.S. Reference Standard Endotoxin (RSE). Endotoxin can be quantified if the concentration is within the range of the standard curve. Water and materials derived from biological sources may contain measurable levels of endotoxin if purification efforts are incomplete. Determined endotoxin content should be compared to the endotoxin limit to assess its significance.

Using the appropriate conditions, Endosafe® Endochrome-K™ has an effective range from 100 - 0.001 EU/mL. Factors influencing the selection of the standard curve range include 1) the parameters of the analytical instrument, 2) the choice of regression models, and 3) the quality of supporting analytical reagents and labware.

BIBLIOGRAPHY

1. Bang, F.B. “A Bacterial Disease of Limulus Polyphemus.” Bull. Johns Hopkins Hosp. 98, p.325 (1956).
2. Cooper, J.F. and Harbert, J.C. “Endotoxins as a Cause of Aseptic Meningitis after Radionuclide Cisternography.” J. Nucl. Med., 16, p.809 (1976).
3. Cooper, J.F., Levin, J., and Wagner, H.N. “Quantitative Comparison of In Vitro and In Vivo Methods for the Detection of Endotoxin.” J. Lab. Clin. Med., 78, p.138 (1971).
4. Hochstein, H.D. “The LAL Test versus the Rabbit Pyrogen Test for Endotoxin Detection: Update ‘87.” Pharm. Technol., 11(6), p.124 (1987).
5. Levin, J. and Bang, F.B. “Clottable protein in Limulus: Its Localization and Kinetics of Its Coagulation by Endotoxin.” Thromb. Diath. Haemorrh., 19, p.186 (1968).
6. McCullough, K.Z. “Process Control: In-process and Raw Material Testing Using LAL.” Pharm. Technol., 12(5) p.40 (1988).
7. Weary, M.E. “Understanding and setting endotoxin limits.” J. Parent. Sci. & Tech., 44:1, p. 16 (1990).
8. Cooper, J.F. “Resolving LAL Test Interferences.” J. Parent. Sci. & Tech., 44:1, p.13 (1990).
9. Guideline on Validation of the Limulus Amebocyte Lysate Test as an End-Product Endotoxin Test for Human and Animal Parenteral Drugs, Biological Products and Medical Devices. U.S. Dept. of Health & Human Services, FDA, December 1987.
10. Interim Guidance for Human and Veterinary Drug Products and Biologicals. U. S. Department of Health & Human Services, FDA, July 15, 1991.
11. “Bacterial Endotoxins Test.” In The U.S. Pharmacopeia, 25<sup>th</sup> revision, 12601 Twinbrook Parkway, Rockville, MD 20852
12. Reference Guide for Endoscan-V, Charles River Laboratories, 1023 Wappoo Road, Suite 43B Charleston, SC, 29407 USA

Manufactured By: CHARLES RIVER ENDOSAFE  
Div. of Charles River Laboratories, Inc.  
1023 Wappoo Road, Suite 43B  
Charleston, SC 29407, USA  
PHONE NUMBER: 843-766-7575  
FAX NUMBER: 843-766-7576

PIR17001