



Chemical Indicators

From indicator tapes and labels to Class 6 emulators, Raven has a full line of chemical indicators for healthcare and industrial users alike. Raven offers indicators for a variety of processes including steam, EO, gamma and dry heat.


ProChem SSI - Steam Sterilization Integrator

- 3 minutes at 270-275°F (132-135°C)
 - Non-reversible color change
- Catalog # CI-SINT
250 Indicators/box
3 year expiry




ProChem EXT - Steam Sterilization Integrator

- Class 5 Integrator
 - 7 minutes at 275°F (135°C)
 - 20 minutes at 250°F (121°C)
 - Extended timing beyond Biological Kill
- Catalog # CI-EXT
250 Indicators/box
5 year expiry



ProChem SSE - Steam Sterilization Emulator

- Class 6 Emulator
 - 3.5 minutes at 275°F (135°C)
 - 4 minutes at 270°F (132°C)
- Catalog # CI-SSIX
250 Indicators/bag
18 month expiry



ProChem OK Cycle 121 - Steam Sterilization Multi-Spot Indicator

- International Use Only*
- Pre-vacuum sterilizers operating at 121°C and 134°C
 - Multiple spots with varied timing
 - Single indicator suitable for multiple cycle lengths
- Catalog # CI-OK121
250 Indicators/bag
3 year expiry



ProChem Process Indicator Tapes

- Tapes for Steam/Chemiclave® and EO
- Distinguishes between processed and unprocessed items
- Good adhesion to woven and non-woven fabrics

Steam

- Catalog # CI-STP-24 24mm x 55M 36 rolls/case
 - Catalog # CI-STP 18mm x 55M 48 rolls/case
 - Catalog # CI-STP-12 12mm x 55M 72 rolls/case
- 3 year expiry

Ethylene Oxide

- Catalog # CI-ETP 18mm x 55M 48 rolls/case
- 2 year expiry



Tape Dispenser
Catalog # CI-T-DISP

ProChem Dry Heat Indicator Labels

- Green to dark brown/black color change
- Distinguishes between processed and unprocessed items



Catalog # CI-DHI
1000 indicator labels/roll
3 year expiry

Gamma Indicator Labels

- Yellow to red color change
- 2 Year Expiry



- Catalog # CI-RO2 - METO® gun label
1000 indicators/roll
- Catalog # CI-RO3 - 12.7mm dots
5000 indicators/roll

Ethylene Oxide Indicator Labels

- Purple to green color change
- 2 Year Expiry



- Catalog # CI-EO2 - METO® gun label
1500 indicators/roll
- Catalog # CI-EO3 - 12.7mm dots
5000 indicators/roll

a division of Mesa Laboratories, Inc.

M E S A



L A B S