



Dry Bath Incubator

Raven Labs offers incubators for all of its biological indicator products. Several configurations are available and incubators are built to suit the needs of your specific area. Incubators maintain a range of +/- 2°C.

- Blocks available to accommodate all Raven BI products (SCBI's¹, ampoules and tubes)
- Built in crusher for SCBI's¹
- Available in 35°C and 55°C²
- Preset temperature eliminates the need for temperature adjustment
- UL listed (120 volt only)
- Plug assemblies available for most countries



ProSpore and ProSpore2 incubators have removable, clear polycarbonate 12-well racks which allows full viewing of all units at a glance. All other incubators operate with aluminum block only. A liquid crystal temperature strip is located on top of the aluminum block (or under the plastic rack for ProSpore and ProSpore2 units) for easy, reliable temperature monitoring.

Incubators available in the following configurations:

Incubators	Number of Wells	120 Volt UL compliant		240 Volt (Outside US)	
		35 Degree	55 ² Degree	35 Degree	55 ² Degree
ProSpore Ampoule 4ml	12 wells	-	I-G55-PS	-	I-240-55-PS
ProSpore Ampoule 1ml	18 wells	-	I-G55-PS1	-	I-240-55-PS1
ProSpore2	12 wells	I-G35-PS2	I-G55-PS2	I-240-35-PS2	I-240-55-PS2
ProTest	14 wells	I-G35-PT	I-G55-PT	I-240-35-PT	I-240-55-PT
ProTest	32 wells	I-G-35-PT32	I-G-55-PT32	I-240-35-PT32	I-240-55-PT32
13mm Test Tubes	12 wells	I-G35-13	I-G55-13	I-240-35-13	I-240-55-13
ProPACT (16mm Test Tubes)	12 wells	I-G35-16	I-G55-16	I-240-35-16	I-240-55-16
ProSpore DH	12 wells	I-G35-DH	-	I-120-35-DH	-
ProAMP	32 wells	-	I-G55-PA	-	I-240-55-PA

¹Self Contained Biological Indicators

²55°C incubators are preset to 58 - 62°C for improved performance of product

