

# Biological Indicators



## PROTEST

- 24-Hour results for Steam; 48-Hour results for EO
- Dramatic color change to or towards yellow
- Easy to handle and activate; fits competitor's incubators

	<u>25 Pack</u>	<u>100 Pack</u>
<b>Steam: Log 5</b>	PT-3-5-25	PT-3-5-100
<b>Log 6</b>	PT-3-6-25	PT-3-6-100
<b>EO: Log 6</b>	PT-1-6-25	PT-1-6-100



## PROSPORE AMPOULE

- For use in Steam liquid loads; suitable for porous loads and hard goods
- Dramatic purple to yellow color change
- 48-Hour results

	<u>4mL ampoule</u>		<u>1mL ampoule</u>	
	<u>10 Pack</u>	<u>50 Pack</u>	<u>15 Pack</u>	<u>100 Pack</u>
<b>Log 4</b>	PS-4-10	PS-4-50	PS1-4-15	PS1-4-100
<b>Log 5</b>	PS-5-10	PS-5-50	PS1-5-15	PS1-5-100
<b>Log 6</b>		PS-6-50	PS1-6-15	PS1-6-100



## PROPACT - BIs for the STERIS System 1®

- Now offering 24-Hour results with new purple to yellow color change
- No refrigeration required; room temperature storage
- Fits competitor's incubators
- Contains 25 test strips, 25 tubes of culture media and 2 transfer clips

Catalog # C-3-PA-24



## PROSPORE2

- 24-Hour results for Steam; 48-Hour results for EO
- Dramatic color change to or towards yellow
- Easy to handle and activate

	<u>50 Pack</u>
<b>Steam: Log 5</b>	PS2-3-5-50
<b>Log 6</b>	PS2-3-6-50
<b>EO: Log 6</b>	PS2-1-6-50



UL compliant incubators available for all self-contained Biological Indicators and test tube sizes (USA models).

Complete catalog available upon request. All products above have received FDA 510k notification and are suitable for sale to the healthcare industry; 510k submission information available upon request. Raven is an FDA registered manufacturer.



# Chemical Indicators

## PROCHEM SSI - Steam Integrator

- Class 5 integrator
- Distinct, non-reversible color change
- 250 strips/box

Catalog # CI-SINT



## PROCHEM EO - Ethylene Oxide Integrator

- Class 5 integrator
- Distinct, non-reversible color change
- 250 strips/box


Catalog # CI-EOINT



## PROCHEM INDICATOR TAPES

- Distinct, non-reversible color change
- Determine immediately if items have or have not been processed

Width	Steam	EO	
24mm	CI-STP-24		(36 Rolls/Case)
18mm	CI-STP	CI-ETP	(48 Rolls/Case)
12mm	CI-STP-12		(72 Rolls/Case)



Dispenser available

## PROCHEM DHI - Dry Heat Indicator Labels

- Distinct, non-reversible color change
- Determine immediately if items have or have not been processed


Catalog # CI-DHI 1000/Roll



## PROCHEM B-D - Air Removal Test Pack

- Lead-free ink
- Clear, non-reversible color change
- Easy to interpret

Catalog # CI-BDTLF20 Box of 20 Tests  
 CI-BDTLF Case of 4 Boxes (80 Tests)



Toll Free: 1-800-728-5702  
 PH: 1-402-593-0781  
 FX: 1-402-593-0995  
 EMAIL: info@ravenlabs.com

**Raven Biological Laboratories**  
 8607 Park Drive  
 Omaha, NE 68127 USA  
[www.ravenlabs.com](http://www.ravenlabs.com)

