

# ProPACT

ProPACT (Peracetic Acid Culture Test kit) is intended for use in monitoring the efficacy of the Steris System 1™\* Processor.



- **24-Hour Results**
- **Ambient storage**
- **Uses existing incubator**
- **Raven incubator available**
- **Contains extra transfer clip**
- **Sample kit available**

## PACKAGING

### Catalog # C-3-PA

#### Contents:

- 25 *G. stearothermophilus* (# 7953<sup>1</sup>) spore strips with a nominal population of  $1.0 \times 10^5$  cfu/strip.
- 25 test tubes (16mm X 75mm) modified Soybean Casein Digest broth recovery media with Bromocresol Purple acid indicator.
- 2 plastic clips to hold strips in place during the cycle.

## CERTIFICATION

Raven's spore strips are certified for purity, resistance and population; Certificate of Analysis accompanies each order. Media is manufactured and tested to meet USP and NCCLS sterility and growth promotion guidelines. ProPACT is stored at 60-80°F, 30-70% RH, away from sterilizing agents, direct sunlight and all other forms of UV light. DO NOT REFRIGERATE. ProPACT has an expiry of one year from date of manufacture. Raven's products are manufactured in accordance with USP Guidelines and the US Food and Drug Administration's Quality System Regulations.

\* Steris and System 1 are registered trademarks of Steris Corporation.

<sup>1</sup> Raven Biological Laboratories derives our spores from reference strains. We do not certify our products as representative of the reference strain and they should therefore not be used where a reference culture is specified.



800.728.5702

Fx 1.402.593.0995

[www.ravenlabs.com](http://www.ravenlabs.com)

[info@ravenlabs.com](mailto:info@ravenlabs.com)

8607 Park Drive

Omaha, NE 68127 USA